

Safecare Biotech (Hangzhou) co.,ltd



SAFECARE COVID-19 Ag

COVID-19 Antigen Rapid Test Device(Swab)



COVID-19 Antigen Rapid Test Device(Swab)

The COVID-19 Antigen Rapid Test is a lateral flow immunoassay intended for qualitative detection of nucleocapsid protein antigen in direct nasal swabs or nasopharyngeal swab specimens from individuals suspected of COVID-19 by their healthcare provider within the first seven days of symptom onset.



SAFECARE COVID-19 Ag Rapid test kit

Intended Use: Detection of SARS-CoV2 Antigen

Package: 25Tests/Box

Storage: 4-30°C

Specimen Type: Nasopharyngeal

Shelf Life: 24 months

Time to Result: 10-15 Minutes

Covid-19 Antigen Rapid Test Device (Swab)

INTRODUCTION




INTENDED USE

For in vitro qualitative detect of Covid-19 antigen in nasal(NS) swab specimens directly from individuals who are suspected of COVID-19 by their healthcare provider within the first 5 days of the onset of the symptoms. This test is only provided for use by clinical laboratories or to healthcare workers for point-of-care testing, and not for at home testing.

PRODUCT PHOTOS





COMPONENT/(BOX)

25 test packed in one box :

25 Test Device

25 Nasal Swab

25 Extraction tube

1 workstation

1 Bottle of buffer

1 Package Insert

PACKAGE SIZE /CARTON

Length:630mm

Width:370mm

Height:300mm

Weight:10Kgs

Include:27 box, 675PCS (product& carton)



SPECIMENS PREPARATION AND TEST PROCEDURE

<Assay Procedure>

Preparation of Extraction solution



Add 6 drops of extraction buffer into the extraction tube



Nasopharyngeal swab collection



Nasal Swab collection



Collect the specimens from the patient using provided swab



Mix the swab with extraction solution



Press the swab and release as much liquid as possible



Place the cap




Add 3 drops of the solution into the sample well



10~15 min



INTERPRETATION OF RESULTS



Positive

Negative

Invalid

CERTIFICATE

**Antigen-Tests zum direkten Erregernachweis des Coronavirus**

safecare

Los

Aktionen

Test-ID	Handelsname	Evaluieru... PEI	Hersteller			Europäischer Bevollmächtigter			Testo...	%
			Name	Stadt	Land	Name	Stadt	Land		
<div style="display: flex; align-items: center; border: 1px solid #ccc; padding: 5px;"> <input checked="" type="checkbox"/> <input style="width: 200px;" type="text" value="Nach 'safecare' suchen"/> </div>										
AT376/21	COVID-19 Antigen Rapid Test Kit (Swab)	Ja	Safecare Biotech (Hangzhou) Co., Ltd.	Hangzhou	CN	NIC GmbH	Osnabrück	DE	POC (ohne Gerät)	
AT483/20	Multi-Respiratory Virus Antigen Test Kit(Swab) (Influenza A+B/ COVID-19)	Ja	Safecare Biotech (Hangzhou) Co., Ltd.	Hangzhou	CN	NIC GmbH	Osnabrück	DE	POC (ohne Gerät)	
AT199/20	COVID-19 Antigen Rapid Test Kit (Swab)	Ja	Safecare Biotech (Hangzhou) Co., Ltd.	Hangzhou	CN	NIC GmbH	Osnabrück	DE	POC (ohne Gerät)	
AT346/21	COVID-19 Antigen Rapid Test Kit (Saliva)	Ja	Safecare Biotech (Hangzhou) Co., Ltd.	Hangzhou	CN	NIC GmbH	Osnabrück	DE	POC (ohne Gerät)	
AT1037/21	COVID-19 Antigen Rapid Test Kit (Swab)	Nein	Safecare Biotech (Hangzhou) Co.,Ltd.	Hangzhou	CN	NIC GmbH	Osnabrück	DE	POC (ohne Gerät)	

CE **CE**

EC Declaration of Conformity

according to the Directive 98/79/EC
(applicable to IVD Devices of NOT Annex II and NOT self-test)

Manufacturer: Safecare Biotech (Hangzhou) Co., Ltd.

Address: Building 2/203, No.18 Haishu Rd.Cangqian Sub-district, Yuhang District, Hangzhou, Zhejiang China 311121

EC Representative: Wellkang Ltd.
16 Castle St,Dover, Kent, CT16 1PW,England,UK
Enterprise Hub,NW Business Complex,1 Beraghmore Road,Derry,BT48 8SE Northern Ireland(NI)

We, the manufacturer, declare under our sole responsibility that

the medical device(s)	Product Name	COVID-19 Antigen Rapid Test Device(Swab)
	Type/model, identification of product allowing traceability (Where applicable)	Cassette(NCO-6012)
of Category:	Common/Others IVD (Devices of NOT Annex II and NOT self-test)	

is/are in conformity with the relevant provisions and requirements of Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

Applied harmonised standards, national standards or other normative documents	EN ISO23640:2015 EN 13612:2002 EN 13641:2002 EN ISO 14971:2019 ISO13485:2016	EN ISO 18113-1:2011 ISO 18113-2: 2009 EN1041- 2008 EN ISO15223-1:2016
---	--	--


Conformity assessment procedure: **Module A (EC Declaration of Conformity) (Annex III, except point e)**

Notified Body (name & number): **NOT applicable**

Certificate & number: Signed on 28th Sep.,2020 Place: Hangzhou, Zhejiang, China

Signature (on behalf of the manufacturer): *Kevin Qiu 2020.9.28*

Name of authorised signatory: Kevin Qiu
Position held in the company: General Manager
Seal Stamp:



Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization
Safecare Biotech (Hangzhou) Co., Ltd.
Building 2/203, No. 18 Haishu Rd,
Cangqian Sub-district, Yuhang District
Hangzhou
311121 Zhejiang
P.R. China

has established and applies a quality management system for medical devices for the following scope:

Design and Development, Manufacture and Distribution of In Vitro Diagnosis of Rapid Test of Fertility, Drug of Abuse, Cardiac Markers, Infectious Diseases


Proof has been furnished that the requirements specified in

EN ISO 13485:2016


are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date:	2020-08-02
Certificate Registration No.:	SX 60149098 0001
An audit was performed. Report No.:	15096152 005
This Certificate is valid until:	2023-06-06

Certification Body



Deutsche
Akkreditierungsstelle
D-91141-01-02



Date 2020-08-02

Herbert Ziegler

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Tel.: +49 221 800-1371 Fax: +49 221 800-3935 e-mail: cert-act@tuev.com http://www.tuev.com/labdy

医疗器械生产许可证

许可证编号:浙食药监械生产许 20140151 号

企业名称: 杭州赛凯生物技术有限公司

生产地址: 杭州市余杭区仓前街道海曙路 18 号 2 号楼
203 室

法定代表人: 唐燕芬

生产范围: 第二类:6840-体外诊断试剂***

企业负责人: 裘科斌

住 所: 杭州市余杭区仓前街道海曙路 18 号 2 号楼
203 室

发证部门: 浙江省药品监督管理局

有效期限: 至 2024 年 8 月 12 日 发证日期: 2019 年 8 月 3 日

1. Easy to collect samples, simple operation, without professional equipment.

2. The test results are available in 15 minutes, and the test results are clearly visible.

3. Convenient transportation and low price, higher accuracy.

4. Suitable for large-scale rapid screening.

ADVANTAGE



SAFECARE BIO-TECH
赛凯生物技术

SAFECARE
COVID-19 Ag



Clinical Evaluation Report

1. Purpose:

In order to verify the clinical performance of the registered test, this clinical evaluation is conducted in R&D lab.

2. Product information:

COVID-19 Antigen Rapid Test Device (Swab) was produced by Safecare Biotech(Hangzhou) Co.,Ltd., Lot number is COV20081201, valid until August,2022.

3. Sample requirement:

Fresh samples were collected from CDC and validated by PCR.

4. Supporting equipment:

PCR tests are performed on ABI7500.
The test-strips are manually operated and visually interpreted.

5. Clinical evaluation:

Researcher: Dr. ZHANG LEI

6. Statistical methods:

		Referencing reagent Test		Total
		Positive	Negative	
Research Reagent	Positive	A	B	A+B
	Negative	C	D	C+D
Total		A+C	B+D	A+B+C+D

Percent Positive Agreement= $A/(A+C)*100\%$

Negative Percent Agreement= $D/(B+D)*100\%$

7. Evaluation indicators:

The total PPA should be no less than 80%.

The total NPA should be no less than 90%.

8. The test data: Refer to the Data Sheet.

9. Statistical results of the clinical evaluation

		Referencing Method (RT-PCR)		Total
		Positive	Negative	
Test-strip	Positive	30	0	30
	Negative	2	52	54
Total		32	52	84

Statistical results

Project	Value	Percentage (95% confidence interval)
Relative Sensitivity (%)	30/32	93.75% (79.19%~99.23%)
Relative Specificity (%)	52/52	100.00% (93.15%~100.00%)

Positive expectation Rate (%)	30/30	100.00% (88.43%~100.00%)
Negative expected Rate (%)	52/54	96.30% (87.25%~99.55%)
Overall Agreement (%)	82/84	97.62% (91.66%~99.71%)

3) Kappa consistency test

According to the literature [5.1], Calculate the Kappa value and standard error; test hypothesis is established for Kappa: $H_0: k = 0$, Kappa value comes from 0 population, $H_1: k > 0$, Kappa value comes from non-0 population, $\alpha = 0.05$.

Project	Value
Kappa Value	0.9489, Good consistency.
Standard Error Se(K)	0.0357
95% Confidence Interval	0.8790~1.0188
Standard Error Se0(K)	0.109
Test Value Z	Z=8.7082, Probability value P=0.0000
Test Result	P<0.05, refuse H_0 , Kappa values come from populations other than 0.

4) Conclusion

A side-by-side comparison was conducted using the research reagent and referencing reagent. Compare with RT-PCR:

The Relative Sensitivity is 93.75%, the Relative Specificity is 100%, the Overall Agreement is 97.62%.

Safecare Biotech (Hangzhou) Co.,Ltd



COMPANY PROFILE

Safecare Biotech(Hangzhou)Co.,Ltd. is a premier and professional manufacturer and supplier of rapid diagnostic test kit with 165 workers, 8000 m² non-dust workshop, a professional R&D team who has 15years experience in rapid test field, advanced automate machines and professional R&D team ensure the high quality, speedy delivery and large production capacity. SAFECARE earned the reputation as a premium brand known for exceptional quality, consistency and innovation.

Our product ranges drug of abuse and alcohol test in urine and saliva, Food Safety test, Women Health test, Infectious Diseases test, Cardiac Markers test and Tumor Markers test with CE & ISO approved. Our drugs tests are even US FDA 510K and CLIA Waived approved which can ensure you high and stable quality.

The available rapid test kits are designed for health-care professionals in laboratories, rehabilitation centers, treatment centers, hospitals, clinics, private practices, human resource departments, mining companies, construction companies and the judicial system. All the products are produced strictly under TUV ISO13485:2016 quality management system for medical devices.

With our highly trained staffs and good service, we are committed to provide professional service and a comprehensive, cutting-edge product offering, help you in selecting the accurate and fast rapid tests and to provide the free samples for your evaluation.

Declaration

Date: 29 Nov. 2021

To whom it may concern,

We, **Safecare Biothech (Hangzhou) Co., Ltd.**, having our office at F1.2 Blog.2, No.18 Haishu Road, Hangzhou 311121 China, hereby declare that the variants listed in the following table can be detected by our **COVID-19 Antigen Rapid Test Kit(Swab) and COVID-19 Antigen Rapid Test Kit(Saliva)**.

Our R&D team has analyzed and designed related recombinant proteins for verification and demonstrate that the products are still effective to the variants.

WHO Label
Alpha
Beta
Gamma
Delta
Omicron

We will keep evaluating the impact of new variants.

Yours sincerely,

杭州赛凯生物技术有限公司
SAFECARE BIOTECH (HANGZHOU) CO., LTD.
HANGZHOU SAFECARE BIOTECH CO.,LTD.