



TECHNICAL DATA SHEET BV-05SC-S Single-use Scrubs

Scrub is fully compliant with EN 13795-1 standards. It is made from a Non-woven SMMS fabric providing AAMI level I-II protection. It offers exceptional protection against fluids under minimal or low liquid pressure conditions.

KEY FEATURES

- •Complies with EN 13795-1:2019.
- •AAMI PB70 level I-II protection.
- •Single-use scrub.
- Latex-free.
- Sterilized with EO
- •Offers protection against biohazards.

STERILE EO

APPROVALS

Fulfill the requirements for Class I MDD under Medical Device Regulation (EU) 2017/745. These scrubs also fulfill all essential requirements and the related rules of the European Economic Commission 93/42/EEC Medical Devices Directives (MDD). They also comply with AAMI Protective Barriers classification system and the associated minimum requirements for the liquid barrier performance of protective appparel and drapes based on ANSI standards.

PACKAGING

Each scrub is individually folded andpackaged into a bag



Veko Giz Plaza 22nd Floor Maslak-İstanbul info@bellevitae.com

Phone: +90 0212 706 1003



SIZE

An appropriate size for the garment should be selected to allow sufficient movement for the wearer whilst maintaining a secure fit



MATERIALS

SUIT

WEIGHT

Non-woven SMMS

43 gr/m2

3XI

67

110

the wearer winist manifalling a secure ht.																							
SIZE CHART																							
MEASURING POINTS (UPPER)		SIZES	1	difference	2	difference	3	difference	4	difference	5	difference	6		MEASURING POINTS		SIZES	1	2	3	4	5	
			S		м		L		XL	unerence	2XL		3XL	Tolerance (+) (-)		(PANTS)	0.2L0	S	Μ	L	XL	2XL	
А	LENGTH	СМ	80	2	82	1.5	83.5	1.5	85	1.5	86.5	1.5	88	1.3	A	WAISTBAND	СМ	57	59	61	63	65	
в	CHEST	СМ	58	2	60	2	62	2	64	2	66	2	68	1.3									
С	ARM LENGTH	СМ	20	1	21	1	22	1	23	1	24	1	25	1.3	в	LENGTH (СМ	105	106	107	108	109	
D	BACK SHOULDER	СМ	51	1	52	1	53	1	54	1	55	1	56	1.3									

INTENDED USE

- •For use in clinical and laboratory settings.
- •Intended to protect from the transfer of microorganisms and bodily fluids in low or minimal risk patient isolation situations.
- •Not intended for high risk surgical procedures, or where the risk of contamination is high.
- •Typical applications may include blood, bodily fluids and environmentally controlled cleanroom

applications.

LIMITED USE





STORAGE AND DISPOSAL

- •Store in dry, clean conditions in original packaging, away from direct sunlight, sources of high temperature, and solvent vapours.
- •Expected shelf life is five years from the date of manufacture when stored as stated.
- •Replace garment if damaged, heavily contaminated or in accordance with local work practice or regulations.
- •Handle and dispose of contaminated garments with care and in accordance with applicable regulations.

WARNING AND LIMITATIONS

Make sure that product is suitable for the application and fitted correctly.Product must never be altered or modified. Final determination as to the suitability of these product for a particular situation is the employer's responsibility. This information is subject to revision at any time.