

TECHNICAL DATA SHEET

BV-07SG Surgical Gown

BV-07SG surgical gown, fully compliant with EN 13795-1 standard, is made from non-woven SMMS fabric providing AAMI level II protection. Level II gowns offer exceptional protection against fluids under minimum or low liquid pressure conditions.

KEY FEATURES



- Complies with EN 13795-1:2019
- AAMI PB70 level II protection
- Single-use gowns
- Ties on the side eliminate the difficulty of tying behind the back
- Velcro closure on the neck provides perfect adjustment for full coverage and flexibility
- Cuffs made of white knitted material on the wrists area provide soft wrist closure for the comfort as well as ensuring a better and tighter fit
- Latex-free
- Non-sterile

NON-STERILE

APPROVALS


Fulfill the requirements for Class I MDD under Medical Device Regulation (EU) 2017/745. These gowns also fulfill all essential requirements and the related rules of the European Economic Commission 93/42/EEC Medical Devices Directives (MDD). They also comply with AAMI Protective Barriers classification system and the associated minimum requirements for the liquid barrier performance of protective apparel and drapes based on ANSI standards



BELLEVITAE
www.bellevitae.com

SIZE

An appropriate size for the garment should be selected to allow sufficient movement for the wearer whilst maintaining a secure fit.



SIZE CHART

SIZE	HEIGHT	CHEST
S	116 cm	65 cm
M	118 cm	68 cm
L	120 cm	71 cm
XL	122 cm	74 cm
XXL	124 cm	77 cm
3XL	126 cm	80 cm

MATERIALS

SUIT Non-woven SMMS

WEIGHT 43 gr/m²

STORAGE AND DISPOSAL

- Store in dry, clean conditions in original packaging, away from direct sunlight, sources of high temperature, and solvent vapours.
- Expected shelf life is five years from the date of manufacture when stored as stated.
- Replace garment if damaged, heavily contaminated or in accordance with local work practice or regulations.
- Handle and dispose of contaminated garments with care and in accordance with applicable regulations.


WARNING AND LIMITATIONS

Make sure that product is suitable for the application and fitted correctly. Product must never be altered or modified. Final determination as to the suitability of these product for a particular situation is the employer's responsibility. This information is subject to revision at any time.

PACKAGING

Each gown is individually folded and packaged into a bag

LIMITED USE



INTENDED USE

- For use in clinical and laboratory settings.
- Intended to protect from the transfer of microorganisms and bodily fluids in low or moderate risk patient isolation situations.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.
- Typical applications may include contact with blood, bodily fluids and environmentally controlled cleanroom applications.

Veko Giz Plaza 22nd Floor Maslak-Istanbul

www.bellevitae.com

info@bellevitae.com

Phone: +90 0212 706 1003