

TECHNICAL DATA SHEET

BV-10SGY Surgical Gown

BV-10SGY Surgical gown is made from non-woven SS fabric providing exceptional protection against fluids under minimum liquid pressure conditions.

KEY FEATURES



- Single-use gowns
- Ties on the side eliminate the difficulty of tying behind the back
- Velcro closure on the neck provides perfect adjustment for full coverage and flexibility
- Cuffs made of white knitted material on the wrists area provide soft wrist closure for the comfort as well as ensuring a better and tighter fit
- Latex-free
- Non-sterile

NON-STERILE

APPROVALS


Fulfill the requirements for Class I MDD under Medical Device Regulation (EU) 2017/745. These gowns also fulfill all essential requirements and the related rules of the European Economic Commission 93/42/EEC Medical Devices Directives (MDD). They also comply with AAMI Protective Barriers classification system and the associated minimum requirements for the liquid barrier performance of protective apparel and drapes based on ANSI standards



BELLEVITAE
www.bellevitae.com

SIZE

An appropriate size for the garment should be selected to allow sufficient movement for the wearer whilst maintaining a secure fit.



SIZE CHART

SIZE	HEIGHT	CHEST
S	116 cm	64 cm
M	119 cm	66 cm
L	122 cm	68 cm
XL	125 cm	71 cm
XXL	128 cm	74 cm
3XL	131 cm	77 cm

MATERIALS

SUIT Non-woven SS

WEIGHT 25 gr/m2

STORAGE AND DISPOSAL

- Store in dry, clean conditions in original packaging, away from direct sunlight, sources of high temperature, and solvent vapours.
- Expected shelf life is five years from the date of manufacture when stored as stated.
- Replace garment if damaged, heavily contaminated or in accordance with local work practice or regulations.
- Handle and dispose of contaminated garments with care and in accordance with applicable regulations.

WARNING AND LIMITATIONS

Make sure that product is suitable for the application and fitted correctly. Product must never be altered or modified. Final determination as to the suitability of these product for a particular situation is the employer's responsibility. This information is subject to revision at any time.

PACKAGING

Each gown is individually folded and packaged into a bag

LIMITED USE



Single use



Do not wash



Do not bleach



Do not iron



Do not dry clean



Do not tumble dry



Flammable

INTENDED USE

- For use in clinical and laboratory settings.
- Intended to protect from the transfer of microorganisms and bodily fluids in minimal risk patient isolation situations.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.