

TECHNICAL DATA SHEET

BV-12SG-R Surgical Gown REINFORCED

BV-12SG-R surgical gown, fully compliant with EN 13795-1 standard, is made from non-woven SMMS fabric with added PE reinforcement providing AAMI level III protection. Level III gowns offer exceptional protection against fluids under moderate or high liquid pressure conditions.

KEY FEATURES



- Complies with EN 13795-1:2019
- AAMI PB70 level III protection
- Single-use gowns
- 2 internal and 2 external double ties on the side eliminate the difficulty of tying behind the back
- Velcro closure on the neck provides perfect adjustment for full coverage and flexibility
- Cuffs made of white knitted material on the wrists area provide soft wrist closure for the comfort as well as ensuring a better and tighter fit
- Anti-static treatment help reduce static build-up.
- PE Reinforcement at the front and forearms
- Latex-free
- Non-Sterile
- Offers protection against biohazards.

NON-STERILE

APPROVALS


Fulfill the requirements for Class I MDD under Medical Device Regulation (EU) 2017/745. These gowns also fulfill all essential requirements and the related rules of the European Economic Commission 93/42/EEC Medical Devices Directives (MDD). They also comply with AAMI Protective Barriers classification system and the associated minimum requirements for the liquid barrier performance of protective apparel and drapes based on ANSI standards



BELLEVITAE
www.bellevitae.com

SIZE

An appropriate size for the garment should be selected to allow sufficient movement for the wearer whilst maintaining a secure fit.



SIZE CHART

SIZE	HEIGHT	CHEST
S	116 cm	65 cm
M	118 cm	68 cm
L	120 cm	71 cm
XL	122 cm	74 cm
XXL	124 cm	77 cm
3XL	126 cm	80 cm

MATERIALS

SUIT	Non-woven SMMS+PE
WEIGHT	55 gr/m ²

STORAGE AND DISPOSAL

- Store in dry, clean conditions in original packaging, away from direct sunlight, sources of high temperature, and solvent vapours.
- Expected shelf life is five years from the date of manufacture when stored as stated.
- Replace garment if damaged, heavily contaminated or in accordance with local work practice or regulations.
- Handle and dispose of contaminated garments with care and in accordance with applicable regulations.

WARNING AND LIMITATIONS

Make sure that the product is suitable for the application and fitted correctly. Product must never be altered or modified. Final determination as to the suitability of these products for a particular situation is the employer's responsibility. Electrostatic Dissipative (ESD) Protective clothing shall not be opened or removed whilst in presence of flammable or explosive atmospheres or while handling flammable or explosive substances. This information is subject to revision at any time.


PACKAGING

Each coverall is individually folded and packaged into a bag.

INTENDED USE

- For use in clinical and surgical operations.
- Intended to protect from the transfer of microorganisms and bodily fluids in moderate or high risk patient surgical situations.
- Typical applications may include contact with blood, bodily fluids and environmentally controlled cleanroom applications.

LIMITED USE



Veko Giz Plaza 22nd Floor Maslak-Istanbul

www.bellevitae.com

info@bellevitae.com

Phone: +90 0212 706 1003