



BELLEVITAE  
www.bellevitae.com



LEVEL 1  
ISOLATION  
GOWN

## TECHNICAL DATASHEET

### BV-13SG Surgical Gown

*BV-13SG Surgical gown is made from non-woven SS fabric providing exceptional protection against fluids under minimum liquid pressure conditions.*

#### KEY FEATURES



- Single-use gowns
- Ties on the side eliminate the difficulty of tying behind the back
- Velcro closure on the neck provides perfect adjustment for full coverage and flexibility
- Cuffs made of white knitted material on the wrists area provide soft wrist closure for the comfort as well as ensuring a better and tighter fit
- Latex-free
- Non-Sterile

#### NON-STERILE

#### APPROVALS

Fulfill the requirements for Class I MDD under Medical Device Regulation (EU) 2017/745. These gowns also fulfill all essential requirements and the related rules of the European Economic Commission 93/42/EEC Medical Devices Directives (MDD). They also comply with AAMI Protective Barriers classification system and the associated minimum requirements for the liquid barrier performance of protective apparel and drapes based on ANSI standards

Veko Giz Plaza 22nd Floor Maslak-İstanbul

[www.bellevitae.com](http://www.bellevitae.com)

[info@bellevitae.com](mailto:info@bellevitae.com)


Phone: +90 0212 706 1003



BELLEVITAE  
www.bellevitae.com

## SIZE

An appropriate size for the garment should be selected to allow sufficient movement for the wearer whilst maintaining a secure fit.



### SIZE CHART

SIZE	HEIGHT	CHEST
S	116 cm	65 cm
M	118 cm	68 cm
L	120 cm	71 cm
XL	122 cm	74 cm
XXL	124 cm	77 cm
3XL	126 cm	80 cm

## MATERIALS

SUIT Non-woven SS

WEIGHT 20 gr/m<sup>2</sup>

## LIMITED USE



Single use



Do not wash



Do not bleach



Do not iron



Do not dry clean



Do not tumble dry



Flammable

## STORAGE AND DISPOSAL

- Store in dry, clean conditions in original packaging, away from direct sunlight, sources of high temperature, and solvent vapours.
- Expected shelf life is five years from the date of manufacture when stored as stated.
- Replace garment if damaged, heavily contaminated or in accordance with local work practice or regulations.
- Handle and dispose of contaminated garments with care and in accordance with applicable regulations.

## WARNING AND LIMITATIONS

Make sure that product is suitable for the application and fitted correctly. Product must never be altered or modified. Final determination as to the suitability of these product for a particular situation is the employer's responsibility. This information is subject to revision at any time.

## PACKAGING

Each gown is individually folded and packaged into a bag

## INTENDED USE

- For use in clinical and laboratory settings.
- Intended to protect from the transfer of microorganisms and bodily fluids in low risk patient isolation situations.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.