

## TECHNICAL DATA SHEET

### BV-17SG Surgical Gown

*BV-17SG Surgical gown is made from non-woven SS fabric providing exceptional protection against fluids under minimum liquid pressure conditions.*

#### KEY FEATURES



- Single-use gowns
- Ties on the side eliminate the difficulty of tying behind the back
- Velcro closure on the neck provides perfect adjustment for full coverage and flexibility
- Cuffs made of white knitted material on the wrists area provide soft wrist closure for the comfort as well as ensuring a better and tighter fit
- Latex-free
- Non-sterile

#### NON-STERILE

#### APPROVALS


Fulfill the requirements for Class I MDD under Medical Device Regulation (EU) 2017/745. These gowns also fulfill all essential requirements and the related rules of the European Economic Commission 93/42/EEC Medical Devices Directives (MDD). They also comply with AAMI Protective Barriers classification system and the associated minimum requirements for the liquid barrier performance of protective apparel and drapes based on ANSI standards



**BELLEVITAE**  
www.bellevitae.com

## SIZE

An appropriate size for the garment should be selected to allow sufficient movement for the wearer whilst maintaining a secure fit.



### SIZE CHART

SIZE	HEIGHT	CHEST
S	116 cm	64 cm
M	119 cm	66 cm
L	122 cm	68 cm
XL	125 cm	71 cm
XXL	128 cm	74 cm
3XL	131 cm	77 cm

## MATERIALS

SUIT	Non-woven SS
WEIGHT	35 gr/m2

## STORAGE AND DISPOSAL

- For use in clinical and laboratory settings.
- Intended to protect from the transfer of microorganisms and bodily fluids in minimal risk patient isolation situations.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.

## WARNING AND LIMITATIONS

Make sure that the product is suitable for the application and fitted correctly. Product must never be altered or modified. Final determination as to the suitability of these products for a particular situation is the employer's responsibility.

## PACKAGING.

Each gown is individually folded and packaged into a bag

## LIMITED USE



Single use



Do not wash



Do not bleach



Do not iron



Do not dry clean



Do not tumble dry



Flammable

## INTENDED USE

- For use in clinical and laboratory settings.
- Intended to protect from the transfer of microorganisms and bodily fluids in minimal risk patient isolation situations.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.

Veko Giz Plaza 22nd Floor Maslak-Istanbul

[www.bellevitae.com](http://www.bellevitae.com)

[info@bellevitae.com](mailto:info@bellevitae.com)

Phone: +90 0212 706 1003